Piezol Implant System

*Surgical Manual
Manuale chirurgico
Chirurgisches Handbuch
Manuel chirurgical
Manual quirúrgico*

Minimally Invasive Technology Simplifying Narrow Ridge Treatment

SURGICAL MANUAL INSTRUCTIONS FOR USE

Caution: Federal law restricts this device to sale only by or on the order of a licensed dentist or physician.

⚠Caution: Federal law restricts this device to sale only by or on the order of a licensed dentist or physician.

General Information and Intended Use

The PiezoImplant System consists of endosseous dental implants, surgical instruments and restorative components in a variety of dimensions to accommodate differing patient anatomy. The REX TL endosseous implants are blade-form having a wedge shape and an apical surface treated with resorbable blast media (grit-blasted and acid passivated). These are offered in buccolingual thicknesses of 1.8mm and 2.9mm having lengths which range from 9mm to 15mm. Cover Screws provide protection to the threads of the abutment connection during the tissue healing process. Retention screws fasten the implant and abutment. A variety of PiezoImplant abutments are offered including Straight, Angled, UCLA, Healing, Cylinder, and Multi-unit. Restorations can be screw and/or cement-retained to the abutments. Laboratory analogs, impression transfers, and pins facilitate creation of the prosthetic restoration. Associated surgical instrumentation includes Piezosurgery® Handpiece inserts for site preparation and the REX IPD with attachments for implantation. A torque of 15 Ncm is recommended for Cover Screws, coping screws, and healing abutments. A torque of 25 Ncm is recommended for retention screws and all other abutments.

The PiezoImplant REX TL implants are provided preassembled to an impression transfer that may also function as a straight abutment via a retention screw to facilitate implantation. Cover Screws are pre-packaged with each preassembled implant. Both the assembled implant and the Cover Screws are provided sterile. The PiezoImplant System surgical instruments and restorative components are provided in a non-sterile condition and must be sterilized prior to use.

The PiezoImplant System implants, abutments, fit gauge, laboratory analogs, impression pins, copings, coping screws, Cover Screws, and retention screws are manufactured from titanium alloy (Ti-6Al-4V ELI) per ASTM F136. UCLA abutments and castable copings are manufactured from polymethylmethacrylate (PMMA). Insert-specific fit gauges are manufactured from stainless steel. The use of restorative components manufactured by Rex Implants is strongly recommended.

The implant platform of the REX TL 1.8 implant series is similar to the 4.1mm diameter (4.1 mmD) external hex implant platform of the Biomet 3i™ Implant System. If abutments from other manufacturers are used with the implants in the REX TL 1.8 series, the angulation must not be greater than 15°, and a Retention Screw TL 1.8 manufactured by Rex Implants must be used to secure the abutment to the implant. The implant platform of the REX TL 2.9 implant series is similar to the 3.5mm diameter (3.5mmD) internal hex implant platform of the Zimmer Dental Screw-Vent® System. If abutments from other manufacturers are used with implants in the REX TL 2.9 series, the angulation must be no greater than 30°.

Indications for Use

The PiezoImplant System is intended for use in dental implant applications for oral rehabilitation of edentulous and partially dentate patients in the maxilla and mandible. Implant retained restorations may consist of single crowns or bridges as well as complete or partial dentures. The prosthetic components are connected to the implants by the corresponding abutments. The PiezoImplant System is intended for delayed loading.

Contraindications

Do not use the Rex PiezoImplant System in presence of following:

- Patients who suffer from any medical conditions that make surgery inadvisable or may be otherwise deemed to be a contraindication by the treating dental practitioner. Such conditions may include and not be limited to: heart disease, diabetes, cirrhosis, HIV positivity, puberty, pregnancy or breastfeeding, radiotherapy, chemotherapy, immunosuppressant therapy, parafunctional therapy, and psychiatric disorder;
- Bone metabolism disturbances,
- Uncontrolled bleeding disorders,
- Inadequate wound healing capacity,
- Uncompleted growth of the maxilla or mandible,
- Drug or alcohol abuse,
- Xerostomia,
- Weakened immune system,
- Uncontrollable endocrine disorders,
- Steroid use,
- Titanium allergy,
- An insufficient quantity of bone volume (height and width) at the implant site,
- Untreated periodontal diseases (loosening of the teeth),
- Untreated severe grinding or clenching of the teeth,
- Infections at the operative site or in the neighboring teeth (pockets, cysts, granulomas), including major sinusitis, or
- Poor hygiene of the mouth and teeth and low compliance (uncooperative or unmotivated).

Relative Contraindications

Caution should be exercised in the presence of the following:

- Exposure to long-term use of opioids or bisphosphonate drugs,
- Previously irradiated bone,
- Diabetes Mellitus,
- Anticoagulation drugs,
- Hemorrhagic diatheses,
- Unfavorable anatomic bone conditions,
- Temporomandibular joint disorders,
- Tobacco use, including moderate to heavy smoking, or
- An unbalanced relationship between the upper and lower teeth.

Warnings and Precautions

- Rex PiezoImplants are single-use products, provided in a triple bag, placed in a carton box and sterilized by gamma irradiation.
- A careful biomechanical study must be performed by the surgeon and restorative dentist to determine the optimal oral restoration for each patient.
- The safe and effective use of Rex PiezoImplants and associated surgical and prosthetic accessories may be achieved only if qualified surgeons trained in the procedure perform the implant treatment.
- A one hundred percent implant success rate cannot be guaranteed.
- The safe and effective use of Rex PiezoImplants, associated surgical instruments, and restorative components is guaranteed only if devices are used as supplied per the instructions provided by Rex Implants and used by properly trained professionals. Modifications to the implants and surgical instruments may result in serious injury or death. Restorative devices may only be modified as directed.
- The surgical and prosthetic procedures required for the appropriate use of these devices are highly specialized and complex. It is essential that the surgical procedure minimizes trauma due to bone overheating and damage to adjacent tissues during the piezoelectric osteotomy.
- Take the following into account when placing a Rex PiezoImplant: 1) the incorrect use or handling of small parts in the patient's oral cavity may lead to inhaling and/or swallowing; and 2) forcing the implant deeper than specifications may lead to expanding or fracturing the implant site bony walls, which could prevent the implant from being secured effectively.
- Techniques required to place and restore dental implants are high complex, requiring specialized knowledge. Practitioners must be trained in implantology techniques before using the PiezoImplant System. Improper surgical technique including a preoperative study that does not properly evaluate the bone volume and the bone quality of the alveolar crest may lead to bone loss and/or implant failure. The clinician should also possess a working knowledge of piezoelectric bone cutting techniques as well as the insertion techniques of press-fit implants.
- Poor surgical planning including failure to reconcile the actual dimensions of surgical instruments to the sizes in radiographic measurements may result in drilling beyond the intended depth. This may result in permanent tissue damage that could result in hemorrhage and/or permanent facial numbness.
- Do not use PiezoImplants in patients if the adequate implant position, sizes, or numbers needed to support functional and eventually parafunctional loading cannot be achieved.
- Primary implant stability is achieved by lightly tapping with the Rex Implant Placement Device (REX IPD) without damaging the bone walls at the surgical site. Primary stability is essential to ensure adequate secondary stability, which is necessary for the long-term success of the implant treatment.
- In the event of implant mobility upon insertion or incorrect implant placement, the implant should be removed. Secure the appropriate Rex Remover Carrier to the implant (Remover Carrier TL 1.8 for the REX TL 1.8 implant series and Remover Carrier TL 2.9 for the REX TL 2.9 implant series). Rex PiezoImplant removal is performed using the REX IPD. The REX IPD is prepared for use by following the instructions in the REX IPD Use and Maintenance Manual. The REX IPD Removal attachment must be secured to the handpiece with the Wrench (do not over-tighten). The Removal action is initiated by selecting the Removal action (↑) on the REX IPD console, setting the force intensity level to 1, and then pressing the Pulse button on the REX IPD console or depressing the REX IPD foot pedal. To remove the implant from the osteotomy, fully engage the Remover Carrier with the REX IPD Remover attachment and initiate the Removal action. During explantation, the REX IPD Remover axis must be parallel with the longitudinal implant axis. If the implant does not move when the Removal action is initiated, incrementally increase the REX IPD force intensity level. Repeat as necessary until the implant is fully removed from the bone. Implant treatment should be repeated only after bone has healed following implant removal. Do not reuse implant and check the Remover Carrier for damage. Do not reuse the Remover Carrier if damage is observed. In cases with a large amount of bone-implant press-fit, the implant may need to be loosened with a piezoelectric cutting device before removal.
- Rex PiezoImplants may be restored only upon completion of the healing process. The REX TL 1.8 system may not be restored until after at least 3 months of healing in high density bone and 6 months in low density bone. The REX TL 2.9 system may not be restored until after at least 4 months of healing in high density bone and 6 months in low density bone.
- Implants placed in the maxilla should not perforate the sinus floor membrane
- During the post-operative healing period, it is critical to protect the implant from trauma and promote osseointegration by ensuring adequate clearance between the implant restoration and antagonist teeth, fixed prosthetic bridge elements or removable prosthetic elements.
- Providing instructions to the patient is essential to ensure success of the implant treatment. The patient must be made aware of the limitations of the implant, importance of oral hygiene and avoidance of contraindications, and potential risk of adverse effects resulting from implant surgery. During the healing period a soft diet must also be prescribed. Patients must be informed to consult a physician if any changes in the performance of the PiezoImplant occurs including bone resorption, loosening, and/or fracture.
- The occlusal load on any Rex PiezoImplant should be similar to that of conventional implants, therefore traumatic and/or parafunctional contacts in the centric relation, right and left laterality, and protrusion must to be avoided.

- The REX TL 1.8 PiezoImplant cannot support a single crown in the molar region even after osseointegration is completed, but can support a single crown if placed elsewhere.
- The REX TL 2.9 PiezoImplant can support a single crown once the osseointegration is completed. If the REX TL 2.9 PiezoImplant is placed in the molar region, a single crown may only be supported in the presence of normal masticatory function.
- If placing more than one REX TL 1.8 or REX TL 2.9 implant in an edentulous area, Rex Implants recommends creating a splint that secures the posts to each other.
- If placing a REX TL 2.9 implant in the mandible, the crestal cortex thickness must be reduced to no more than 2mm.
- If abutments from other manufacturers are used with the implants in the REX TL 1.8 series, the angulation must not be greater than 15°, and a Retention Screw TL 1.8 manufactured by Rex Implants must be used to secure the abutment to the implant.
- If abutments from other manufacturers are used with implants in the REX TL 2.9 series, the angulation must be no greater than 30°.
- After implant placement, if more than 40% of a REX TL PiezoImplant surface is exposed (i.e. not surrounded by bone), the implant should be removed from the patient.

Potential adverse events

Potential adverse events must be communicated to the patient prior to surgery. Potential adverse events related to the use of dental implants may include: integration failure; integration reduction; wound dehiscence requiring bone graft; jaw bone fracture; the perforation of the following: maxillary sinus, inferior border of the mandible, labial and lingual bone walls, alveolar canal, and gingiva; abscesses, fistulas, suppuration, inflammation, radiotransparency, persistent pain, sensitivity reduction, paresthesia, hyperplasia, excessive bone reduction requiring surgery, implant fracture, systemic infections, nerve lesions or other nerve damage, and vascular lesions or hemorrhaging, which at times may be serious especially in patients undergoing treatment with anticoagulants and/or antiaggregants.

Failure to comply with the instructions contained in this document, including reuse of products marked as single-use, may cause harm to the patient, including risk of serious injury or death.

How supplied

The REX TL PiezoImplants are provided preassembled to a transfer, which may also function as a straight abutment, via a retention screw to facilitate implantation. Cover Screws are pre-packaged with each preassembled implant.

The implant assembly and included Cover Screw are provided sterile and should be handled in a safe manner. All sterile products are labeled “STERILE.” Do not use the sterile products if the package is damaged or already opened. Do not use the sterile products if the expiration date on the label has passed. Do not clean, re-sterilize or reuse any sterile product. Due to the complex surface characteristics of the PiezoImplants, only gamma radiation has been demonstrated to effectively sterilize these devices. The REX TL PiezoImplants should be stored in a clean, dry, dust-free, dark place at room temperature.

The PiezoImplant System restorative components, surgical tray, and surgical instruments are provided non-sterile and must be sterilized prior to use per the instructions provided below. Piezosurgery® inserts, the REX IPD, and REX IPD accessories are also provided non-sterile and must be cleaned and sterilized per the instructions for use provided with the device. All dental abutments and restorative devices are indicated as single use.

Cleaning/Reprocessing

The devices in the PiezoImplant System provided non-sterile must be thoroughly cleaned before initial use. Reusable devices must be cleaned between uses. While it is recommended that the following validated steps are included in a reprocessing protocol, the end-user bears the ultimate responsibility for the cleanliness of the device. These instructions are not intended for implants or devices not manufactured by Rex Implants.

Manual Cleaning/Disinfecting Reprocessing Method

1. Rinse device in cold, potable tap water (< 43°C; <109°F) to remove debris and prevent coagulation of blood.
2. Prepare a solution of enzyme detergent and potable tap water at pH 7, according to the manufacturer's instructions.
3. Place the device in a container. Add a sufficient amount of solution of enzyme detergent in the container to cover completely.
4. Leave the device immersed for 10 minutes at 40°C (104°F). This will reduce the organic residues. When immersed in the solution of enzyme detergent, gently brush any threaded parts and small grooves on the device.
5. Place the device in an ultrasonic bath containing enzymatic detergent solution at 40°C (104°F) for at least 10 minutes. This will reduce the organic compounds on the device. After soaking in the ultrasonic bath, gently brush any threaded parts or grooves on the device.
6. Under running, warm, potable tap water and using a toothbrush with soft bristles made of nylon, clean the device thoroughly without damaging the surface. Do this until all visible traces of soil are removed. Carry out the final rinse with distilled water.
7. Finish cleaning and inspect the devices under an appropriate light source, paying attention to details that might house soil (e.g., threads, holes, slots) and if necessary repeat the cleaning cycle.
8. Inspect all devices for any signs of wear and tear. Do not use any device whose integrity is visibly compromised.
9. Dry the device in preparation for sterilization.

Automated Cleaning/Disinfecting Reprocessing Method

1. Lay the device in a metallic tray and place it in thermodisinfectior.
- NOTE: Place the devices in the washing machine so that dead zones do not arise and the water can properly drain. Also, make sure

that the devices are properly held in place in the washing basket and cannot move during the washing process, as shocks could damage them.

WARNING: Avoid overloading. Overloading compromises cleaning effectiveness.

2. Set the following sequence and parameters applicable to the cycle:

- 3 min, Pre-wash with cold, potable tap water;
- 5 min, Wash with enzymatic cleaner at $45^{\circ}\text{C} \pm 2^{\circ}\text{C}$ (113°F);
- 2 min, Rinse with cold, potable tap water;
- 5 min, Thermodisinfest at 93°C (200°F) with demineralized water following the national requirements about A0 values (thermal disinfection at a temperature of 90°C [200°F] for 5min results in an A0 value of 3000);
- 2 min rinse with cold, de-mineralized water.

3. Finish cleaning and inspect the devices under an appropriate light source, paying attention to details that might house soil (e.g., threads, holes, slots) and if necessary repeat the cleaning cycle.

4. Inspect all devices for any signs of wear and tear. Do not use any device whose integrity is visibly compromised.

5. Dry the device in preparation for sterilization.

Drying

Before starting the sterilization cycle, make sure that the device is thoroughly dry, both externally and internally. For this purpose, blow compressed air both externally and into/through any holes; this will prevent the onset of stains, haloes, or rust on the device.

Sterilization

The PiezoImplant System abutments, surgical tray, surgical instruments, and restorative accessories are provided non-sterile and must be cleaned and sterilized prior to use. These instructions are not intended for implants or devices not manufactured by Rex Implants such as Piezosurgery® Inserts. Prior to sterilization, wrap all instruments individually using FDA cleared sterilization pouches. Drying is to be performed inside the steam sterilizer in the prevacuum cycle at 132°C (270°F). All sterilization phases must be performed by the operator in compliance to ANSI/AAMI/ISO 17665-1:2006, EN ISO 556-1:2001, and ANSI/AAMI ST99:2010.

Sterilization validation has shown the following recommendations for sterilization to be effective to an SAL of 10^{-6} .

Method: Steam

Cycle: Prevacuum

Temperature: 132°C (270°F)

Exposure Time: 4 minutes

Minimum Drying Time: 20 minutes

MRI Safety Information

Non-clinical testing has demonstrated the Rex PiezoImplant is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 T and 3 T,
- Maximum spatial field gradient of 1700 gauss/cm, and
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode). Under the scan conditions defined above, the Rex PiezoImplant is expected to produce a maximum temperature rise of less than 2.2°C after 15 minutes of continuous scanning. In non-clinical testing, the image artifact caused by the device extends approximately 16.2 mm from the Rex PiezoImplant when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

Opening the Sterile Packaging

Open the cardboard box and extract its contents using non-sterile gloves. Manually open the Mylar bag in the direction of the Chevron peel opening side and drop the sterile Tyvek pouch into a sterile tray. Do not allow the Tyvek pouch to contact any non-sterile objects. Use sterile gloves to open the Tyvek pouch in the direction of the Chevron peel opening side and pull out the polybag pouch. Hold the packaged, pre-assembled implant and open the larger portion of the pouch using surgical scissors without dropping the contents. Connect the Thumb Knob (TK) to the square head of the transfer that is pre-assembled to the implant and extract the pre-assembled implant from the polypouch. After the implant has been placed, hold the packaged Cover Screw and open the smaller portion of the pouch using surgical scissors without dropping the contents. Remove the Cover Screw from the pouch.

Surgical Manual

The proper implant placement procedure consists of several surgical steps, which are briefly described below. The instructions provided below are general guidelines and may not be applicable to every surgical case. The surgeon must determine the actual steps comprising the surgical procedure by applying his or her best clinical judgment to the study of each patient's unique anatomy.

Step 1: Pre-operative Planning

Techniques required to place and restore dental implants are highly complex, requiring specialized knowledge. A careful biomechanical study must be performed by the surgeon and restorative dentist to determine the optimal oral restoration for each patient. Practitioners must be trained in implantology techniques before using the PiezoImplant System to achieve safe and effective use of these devices. Improper surgical technique, including a preoperative study that does not properly evaluate the bone volume and the bone quality of the alveolar crest, may lead to bone loss and/or implant failure. The clinician should also possess a working knowledge of piezoelectric bone cutting techniques as well as the insertion techniques of press-fit implants.

A PiezoImplant radiographic template for the pre-surgical study of each implant site is provided. The implant radiographic template should be used in combination with cross section radiographic images or Cone Beam CT scans in order to evaluate the thickness and quality of the residual crestal bone, to evaluate the proper position of the implant site, and to determine the most suitable PiezoImplant size to be used (Figure 1). Do not place implants in the presence of local root remains.

Do not use PiezoImplants in patients if the adequate implant position, sizes, or numbers needed to support functional and eventually parafunctional loading cannot be achieved. Poor surgical planning, including failure to reconcile the actual dimensions of surgical instruments to the sizes in radiographic measurements, may result in preparing the site beyond the intended depth. This may result in permanent tissue damage that could result in hemorrhage and/or permanent facial numbness.

PiezoImplants are recommended for the treatment of thin residual crestal bone at least 9 - 16 mm high and 3.5 mm wide for the REX TL 1.8 implant series and 5 mm wide for the REX TL 2.9 implant series. A distance of 1 mm is required between the REX TL PiezoImplants and other implants or proximal teeth. Typical osteotomies for various bone conditions are available in Tables 1-4 for the REX TL 1.8 implant series and Tables 5-8 for the REX TL 2.9 implant series. These tables do not provide clinical recommendations, and the surgeon must use his or her expert opinion to determine the best osteotomy preparation for each patient on an individual basis.




Figure 1. Pre-operative planning (REX TL 1.8 Radiographic Template used)

REX TL 1.8

Table 1. Typical osteotomy descriptions for **REX TL 1.8 H9** implants placed in various bone conditions

Implant	Bone Type	Extremely Hard Crestal Bone Depth (mm)	Osteotomy Description (Laser Marking Reached)				
			W1	W2	W3	W4	W4-H
REX TL 1.8 H9	Normal/Hard	0	2	1	1	1	N/A
		1	2	1	1	1	1
	Hard	2	2	1	1	1	2
		3	2	1	1	1	3

Table 2. Typical osteotomy descriptions for **REX TL 1.8 H11** implants placed in various bone conditions

Implant	Bone Type	Extremely Hard Crestal Bone Depth (mm)	Osteotomy Description (Laser Marking Reached)				
			W1	W2	W3	W4	W4-H
REX TL 1.8 H11	Normal	0	2	2	2	1	N/A
		1	2	2	2	1	1
		2	2	2	2	1	2
		3	2	2	2	1	3
	Hard	0	2	2	2	2	N/A
		1	2	2	2	2	1
		2	2	2	2	2	2
		3	2	2	2	2	3

Table 3. Typical osteotomy descriptions for **REX TL 1.8 H13** implants placed in various bone conditions

Implant	Bone Type	Extremely Hard Crestal Bone Depth (mm)	Osteotomy Description (Laser Marking Reached)				
			W1	W2	W3	W4	W4-H
REX TL 1.8 H13	Normal	0	2	3	3	1	N/A
		1	2	3	3	1	1
		2	2	3	3	1	2
		3	2	3	3	1	3
	Hard	0	2	3	3	2	N/A
		1	2	3	3	2	1
		2	2	3	3	2	2
		3	2	3	3	2	3

Table 4. Typical osteotomy descriptions for **REX TL 1.8 H15** implants placed in various bone conditions

Implant	Bone Type	Extremely Hard Crestal Bone Depth (mm)	Osteotomy Description (Laser Marking Reached)				
			W1	W2	W3	W4	W4-H
REX TL 1.8 H15	Normal	0	2	4	4	1	N/A
		1	2	4	4	1	1
		2	2	4	4	1	2
		3	2	4	4	1	3
	Hard	0	2	4	4	2	N/A
		1	2	4	4	2	1
		2	2	4	4	2	2
		3	2	4	4	2	3

REX TL 2.9

Table 5. Typical osteotomy descriptions for **REX TL 2.9 H9** implants placed in various bone conditions

Implant	Bone Type	Extremely Hard Crestal Bone Depth (mm)	Osteotomy Description (Laser Marking Reached)							
			W1	W2	W3	W4	W4-H	W5	W6	W6-H
REX TL 2.9 H9	Normal	0	2	1	1	1	N/A	1	N/A	N/A
		1	2	1	1	1	2	1	N/A	1
		2	2	1	1	1	2	1	1	2
		3	2	1	1	1	3	1	2	3
	Hard	0	2	1	1	1	N/A	1	2	N/A
		1	2	1	1	1	3	1	2	1
		2	2	1	1	1	3	1	2	2
		3	2	1	1	1	3	1	2	3

Table 6. Typical osteotomy descriptions for **REX TL 2.9 H11** implants placed in various bone conditions

Implant	Bone Type	Extremely Hard Crestal Bone Depth (mm)	Osteotomy Description (Laser Marking Reached)							
			W1	W2	W3	W4	W4-H	W5	W6	W6-H
REX TL 2.9 H11	Normal	0	2	2	2	2	N/A	2	N/A	N/A
		1	2	2	2	2	3	2	N/A	1
		2	2	2	2	2	3	2	1	2
		3	2	2	2	2	3	2	2	3
	Hard	0	2	2	2	2	N/A	2	2	N/A
		1	2	2	2	2	3	2	2	1
		2	2	2	2	2	3	2	2	2
		3	2	2	2	2	3	2	2	3


Table 7. Typical osteotomy descriptions for **REX TL 2.9 H13** implants placed in various bone conditions

Implant	Bone Type	Extremely Hard Crestal Bone Depth (mm)	Osteotomy Description (Laser Marking Reached)							
			W1	W2	W3	W4	W4-H	W5	W6	W6-H
REX TL 2.9 H13	Normal	0	2	3	3	3	N/A	3	N/A	N/A
		1	2	3	3	3	3	3	N/A	1
		2	2	3	3	3	3	3	1	2
		3	2	3	3	3	3	3	2	3
	Hard	0	2	3	3	3	N/A	3	3	N/A
		1	2	3	3	3	3	3	3	1
		2	2	3	3	3	3	3	3	2
		3	2	3	3	3	3	3	3	3

Table 8. Typical osteotomy descriptions for **REX TL 2.9 H15** implants placed in various bone conditions

Implant	Bone Type	Extremely Hard Crestal Bone Depth (mm)	Osteotomy Description (Laser Marking Reached)							
			W1	W2	W3	W4	W4-H	W5	W6	W6-H
REX TL 2.9 H15	Normal	0	2	4	4	4	N/A	4	N/A	N/A
		1	2	4	4	4	3	4	N/A	1
		2	2	4	4	4	3	4	1	2
		3	2	4	4	4	3	4	2	3
	Hard	0	2	4	4	4	N/A	4	4	N/A
		1	2	4	4	4	3	4	4	1
		2	2	4	4	4	3	4	4	2
		3	2	4	4	4	3	4	4	3

Step 2: Resect soft tissue and make proper measurements

Resect the soft tissue as necessary to expose the implantation site. Measure the available ridge thickness for implant placement with a standard dental probe. The ridge width should be no less than 3.5 mm and between 9 and 16 mm high. In the case of partial edentulism, the pilot osteotomy should be created 8 mm from the axis of the last single-rooted tooth (Figure 1) or other dental implant. In the case of mono-edentulism, the pilot osteotomy should be created 8 mm mesio-distally from the axis of proximal teeth or 1.5 mm distal of the root to maintain a minimum of 1.5 mm between implant and root.

Figure 1: Resect soft tissue and make measurements

Step 3: Create pilot osteotomy

Osteotomies are only to be created using a Piezosurgery® device by Mectron. Prepare a Piezosurgery® device for use by following the instructions for use provided with the device. A cut is complete when the insert moves freely within the osteotomy.

Attach the W1 Piezosurgery® insert to the Piezosurgery® handpiece. Create the initial pilot osteotomy at the center of the desired implant position, leaving a space of at least 1.5 mm between the implant and any remaining teeth (Figure 2). The correct depth of osteotomy (9 mm) is ensured when the bottom of the larger, upper laser marking is reached.


Figure 2: Create pilot osteotomy (W1 insert)

Step 4: Verify pilot osteotomy axis

Use the alignment pin to verify the proper position and angulation of the osteotomy and check the distance from the proximal teeth.



Figure 3: Verify axis of pilot osteotomy (Alignment Pin)

Step 5: Create reference pilot osteotomies

Use the W1 insert to create a reference pilot osteotomy on each side (mesial and distal) of the initial pilot osteotomy at a distance of 1 mm from the initial osteotomy. The osteotomy is considered complete upon reaching the bottom of the lower, smaller laser marking (Figure 3). The total distance between the reference osteotomies must be at least 5 mm and not exceed 6.5 mm.


Figure 4: Create reference pilot osteotomy (W1 insert)

Step 6: Connect pilot osteotomies

Attach the W2 micro-saw insert tip to the Piezosurgery® handpiece. Connect the pilot osteotomies using the W2 insert (Figure 5). The depth of the osteotomy to be created (9, 11, 13, 15mm) depends on the length of the implant, reflected in the laser markings on the insert. Laser markings are provided to designate depths of 9 mm (bottom/apical marking), 11, 13, and 15 mm (top/coronal marking). The depth of the osteotomy should not exceed the length of the PiezoImplant to be used.


Figure 5: Connect pilot osteotomies (W2 insert)

Step 7: Refine the osteotomy

Attach the W3 micro-file insert tip to the Piezosurgery® handpiece. Refine the osteotomy with the W3 insert (Figure 6), stopping when the laser markings corresponding to the length (9, 11, 13, 15mm) of the PiezoImplant to be used is reached.


Figure 6: Refine the osteotomy (W3 insert)

Step 8: Verify proper osteotomy length and alignment

Verify the dimensions and alignment of the osteotomy with the Fit Gauge, which must be inserted to the laser mark corresponding to the length of the Rex PiezoImplant to be used (Figure 7). The Fit Gauge has laser marks indicating 9 mm (bottom/apical mark), 11 mm, 13 mm, and 15 mm (top/coronal mark), corresponding to the lengths of the REX TL implants. The appropriate laser markings should be slightly below bone level. If the Fit Gauge does not fit within the osteotomy, repeat steps 6 and 7.


Figure 7: Verify proper osteotomy length and alignment (Fit Gauge)

Step 9: Enlarge the osteotomy buccolingually

Attach the W4 wedge-shaped, micro-file insert tip to the Piesosurgery® handpiece. Enlarge the osteotomy buccolingually in the crestal bone using the W4 insert (Figure 8). The W4 insert has a maximum width of 1.6 mm and has laser markings at 3, 5, 7, and 9 mm of working length for reference. The correct depth must be determined by the surgeon after a careful study of the patient anatomy including bone density.



Figure 8: Enlarge the osteotomy mesiodistally buccolingually (W4 insert shown)

During the osteotomy, frequently verify the proper degree of underpreparation with the Fit Gauge W4/W4-H, by inserting its long, stepped end in the osteotomy (labeled W4). Laser markings indicate preparation depths of 3, 5, 7, and 9 mm.

WARNING: Failure to provide proper buccolingual space in the osteotomy increases the risk of iatrogenic bone fracture.

REX TL 1.8 series implants: The depth of the W4 Insert should not exceed the top of the laser mark at the base of the wedge (3 mm) for most cases. For hard, mineralized bone, the depth of the osteotomy may be increased by continuing the osteotomy. For implants 9 mm in length, do not exceed the first laser marking (3 mm). For implants 11 mm in length or greater, do not exceed the second laser marking (5 mm).

REX TL 2.9 series implants: Enlarge the osteotomy with the W4 insert to the desired osteotomy depth as indicated by the laser markings corresponding to the length (9, 11, 13, 15 mm) of the PiezoImplant to be used. For implants 9 mm in length, do not exceed the first laser marking (3 mm).

Step 10: (Optional) – In extremely hard bone, continue to enlarge the osteotomy buccolingually in the coronal region

For extremely hard, mineralized bone, attach the W4-H micro-file insert tip to the Piesosurgery® handpiece. Enlarge the osteotomy buccolingually in the crestal bone using the W4-H insert (Figure 9). The W4-H insert is 1.8 mm wide and has laser markings at 1, 2, and 3 mm of working length for reference. The correct osteotomy depth must be determined by the surgeon after a careful study of the patient anatomy. Never exceed the third laser marking (3 mm) for any implant length.



Figure 9: For extremely hard, mineralized bone continue to enlarge the osteotomy buccolingually (optional; W4-H insert shown)

During the osteotomy, frequently verify the proper degree of underpreparation with the Fit Gauge W4/W4-H by inserting its short, straight end (labeled W4-H) in the osteotomy. Laser markings indicate preparation depths of 1, 2, and 3 mm

WARNING: Use of the W4-H insert creates an osteotomy that is the same size as the coronal, intraosseous portion of the implant, thus removing compression of the implant within the bone. Exceeding the third laser marking of the insert may result in lack of primary stability, making immediate REX TL 1.8 implant placement impossible.

Step 11: REX TL 2.9 SERIES ONLY – Continue to enlarge the osteotomy

If a REX TL 2.9 series PiezoImplant is to be placed in the osteotomy attach the W5, micro-file insert tip to the Piezosurgery® handpiece. Continue to enlarge the osteotomy to the desired size using the W5 insert as indicated by the laser markings on the shaft corresponding to the length (9, 11, 13, 15 mm) of the PiezoImplant to be used. The W5 insert is 2.3 mm wide and has laser markings at 2, 3, 4, and 5 mm of working length for reference.

During the osteotomy, frequently verify the proper degree of underpreparation with the Fit Gauge W5. Laser markings indicate preparation depths of 2, 3, 4, and 5 mm.

For hard, mineralized bone, continue to enlarge the osteotomy with the W6 insert as determined to be necessary by pre-surgical planning. The W6 insert is 2.7mm wide and has laser markings at 2, 3, 4, and 5 mm of working length for reference. For implants 9 and 11 mm in length, do not exceed the second laser marking (3 mm). For implants 13 mm in length, do not exceed the third laser marking (4 mm).

Never exceed the fourth laser marking (5 mm) for any implant length.

During the osteotomy, frequently verify the proper degree of underpreparation with the Fit Gauge W6/W6-H, by inserting its longer end (labeled W6) in the osteotomy. Laser markings indicate preparation depths of 2, 3, 4, and 5 mm.

Figure 10: For the REX TL 2.9 series only, continue to enlarge the osteotomy mesiodistally buccolingually (W5 insert shown)

Step 12: (Optional) for REX TL 2.9 series only – In extremely hard bone, continue to enlarge the osteotomy buccolingually in the coronal region

For extremely hard, mineralized bone, attach the W6-H micro-file insert tip to the Piezosurgery® handpiece. Enlarge the osteotomy buccolingually in the crestal bone using the W6-H insert (Figure 11). The W6-H insert is 2.9 mm wide and has laser markings at 1, 2, and 3 mm of working length for reference. The correct osteotomy depth must be determined by the surgeon after a careful study of the patient anatomy. Never exceed the third laser marking (3mm) for any implant length.

During the osteotomy, the surgeon must frequently verify the proper degree of underpreparation with the Fit Gauge W6/W6-H, by inserting its short end (labeled W6-H) in the osteotomy. The markings on this end of the Fit Gauge W6/W6-H indicate 1, 2, and 3 mm of preparation depth.

WARNING: Use of the W6-H insert creates an osteotomy that is the same size as the coronal, intraosseous portion of the implant, thus removing compression of the implant within the bone. Exceeding the third laser marking of the insert may result in lack of primary stability and make placing a REX TL 2.9 implant impossible.


Figure 11: For extremely hard, mineralized bone continue to enlarge the osteotomy buccolingually (optional; W6-H insert shown)

Step 13: Remove the PiezoImplant for the sterile packaging and place the PiezoImplant it in the osteotomy site.

Use sterile gloves to open the Tyvek implant pouch in the direction of the Chevron peel opening side and pull out the polybag pouch. Hold the packaged, pre-assembled implant and open the larger portion of the pouch using surgical scissors without dropping the contents. Engage the self-retaining clips of the Thumb Knob over the square head of the pre-assembled transfer/straight abutment (Figure 12, top), and carefully withdraw the assembly from the packaging.

Place the implant into the osteotomy site by applying a small force to the Thumb Knob (Figure 12, bottom). Detach the Thumb Knob from the PiezoImplant Assembly.

Note: The Cover Screw is packaged in a separate pouch compartment with the PiezoImplant assembly. Ensure that the Cover Screw pouch is not inadvertently opened when opening the PiezoImplant pouch and maintain the packaging sterility.




Figure 12: Extract the PiezoImplant assembly from packaging (top; REX TL 1.8 implant series shown). Place the PiezoImplant into the osteotomy site (bottom).

Step 14: Use the REX IPD to seat the implant in the osteotomy

Use the Rex Implant Placement Device (REX IPD) to fully insert the PiezoImplant into the osteotomy site (Figure 13). The following description is not sufficient to allow immediate use of the REX IPD. Please refer to the REX IPD Use and Maintenance Manual for complete instructions. The REX IPD is prepared for use by following the instructions in the REX IPD Use and Maintenance Manual.

Connect the keyboard console, handpiece, and optional foot pedal if desired. Turn the power switch on. Select the appropriate Striker attachment (anterior or posterior) and secure to the REX IPD handpiece with the Wrench (do not over-tighten). The intensity of the insertion action is controlled by the intensity buttons (1 through 4). The insertion action is initiated by selecting the Strike action (↓) and lowest force intensity value on the REX IPD console, and then pressing the Pulse button on the REX IPD console or depressing the REX IPD foot pedal.

Seat the striker tip into the head of the Transfer/Straight Abutment. Press the Pulse button (or press the REX IPD foot pedal) to impact the implant. If the implant does not move when the foot pedal / Pulse button is pressed, incrementally increase the intensity level on the REX IPD console. Repeat as necessary until the implant is fully seated with the shoulder of the implant slightly beneath the bone level as shown in Figure 14. When the implant is fully seated into the osteotomy, disengage the REX IPD striker.

In the event of implant mobility upon insertion or incorrect implant placement, the implant should be removed immediately by first securing the appropriate Rex Removal Carrier to the implant (Remover Carrier TL 1.8 for the REX TL 1.8 implant series and Remover Carrier TL 2.9 for the REX TL 2.9 implant series).

Rex PiezoImplant removal is performed using the REX IPD. The REX IPD Removal attachment must be secured to the handpiece with the Wrench (do not over-tighten). The Removal action is initiated by selecting the Removal action (↑) on the IPD console, setting the force intensity level to 1, and then pressing the Pulse button on the IPD console or depressing the IPD foot pedal. To remove the implant from the osteotomy, fully engage the Remover Carrier with the REX IPD Remover attachment and initiate the Removal action.

During explantation, the REX IPD remover axis must be parallel with the longitudinal implant axis. If the implant does not move when the Removal action is initiated, incrementally increase the REX IPD force intensity level. Repeat as necessary until the implant is fully removed from the bone. Implant treatment should be repeated only after bone has healed following implant removal. Do not reuse implant.

WARNING: Do not use the REX IPD and/or Remover Carrier for implant removal after healing has begun.



Figure 13: Use the REX IPD to seat the implant in the osteotomy (IPD Striker, Assembled PiezoImplant – REX TL 1.8 series shown).

Step 12: Remove the Transfer / Straight Abutment

To remove the Transfer / Straight Abutment, use a 0.050" hex driver with the thumb knob (TK) and use the assembly to disengage the Retention Screw.


Figure 12: Remove the Transfer/Straight Abutment from the PiezoImplant (REX TL 1.8 series shown) with a driver

Step 16: Attach a Cover Screw and Suture Soft Tissue

Hold the packaged Cover Screw and open the smaller portion of the poly-bag pouch using surgical scissors without dropping the contents. Remove the Cover Screw from the pouch and hand-tighten to the implant. Connect a 0.050" hex driver to the Cover Screw and tighten with 15 Ncm of torque. Remove the driver from the Cover Screw (Figure 15). Additional implants are placed as necessary by repeating the previous instructions. The soft tissue is sutured around the Cover Screw.

Additional implants are placed as necessary by repeating the previous instructions. The soft tissue is sutured around the Cover Screw.


Figure 15: Attach the Cover Screw to the PiezoImplant (REX TL 1.8 series shown) with a torque of 15 Ncm

Step 17: Post-operative Healing Period

During the post-operative healing period, it is critical to protect the implant from trauma and promote osseointegration by ensuring adequate clearance between the implant restoration and antagonist teeth, fixed prosthetic bridge elements, or removable prosthetic elements. If placing more than one REX TL 1.8 or REX TL 2.9 implant in an edentulous area, Rex Implants recommends creating a splint that secures the posts to each other.

A PiezoImplant System Patient Passport is provided with each implant. Place a patient label for each implant used in the appropriate location on this document. The Patient Passport serves as a patient record and must be given to the patient or designated caregiver after implant placement surgery is completed.

Providing instructions to the patient is essential for successful implant treatment. All patients must be made aware of the limitations of the PiezoImplant System, importance of oral hygiene and avoidance of contraindications, and potential risk of adverse effects resulting from implant surgery. During the healing period a soft diet must also be prescribed. Patients must be informed to consult a physician if any changes in the performance of the PiezoImplant occurs including bone resorption, loosening, and/or fracture.

The REX TL 1.8 system may not be restored until after at least 3 months of healing in high density bone and 6 months in low density bone. The REX TL 2.9 system may not be restored until after at least 4 months of healing in high density bone and 6 months in low density bone.

Step 18: Application of Prosthetic Components

Rex PiezoImplants may be restored only upon completion of the healing process. Once osseointegration of the implant has been achieved, the prosthetic component(s) are applied in a secondary procedure. A careful biomechanical study should be performed to determine the optimal restoration for each patient, and though PiezoImplants have a history of clinical success, a one hundred percent success rate cannot be guaranteed.


The REX TL 1.8 PiezoImplant series cannot support a single crown in the molar region even after osseointegration is completed, but can support a single crown if placed elsewhere. The REX TL 2.9 PiezoImplant series can support a single crown once the osseointegration is completed. If the REX TL 2.9 PiezoImplant is placed in the molar region, a single crown may only be supported in the presence of normal masticatory function.

The implant platform of the REX TL 1.8 implant series is similar to the 4.1mm diameter (4.1 mmD) external hex implant platform of the Biomet 3i® Implant System. If abutments from other manufacturers are used with the implants in the REX TL 1.8 series, the angulation must not be greater than 15°, and a Retention Screw TL 1.8 manufactured by Rex Implants must be used to secure the abutment to the implant. The implant platform of the REX TL 2.9 implant series is similar to the 3.5mm diameter (3.5mmD) internal hex implant platform of the Zimmer Dental Tapered Screw-Vent® System. If abutments from other manufacturers are used with implants in the REX TL 2.9 series, the angulation must be no greater than 30°.

Though restorative components from other manufacturers may fit on the REX TL 1.8 PiezoImplant series and REX TL 2.9 PiezoImplant series, the use of restorative components manufactured by Rex Implants is strongly recommended. Instruct the patient to bring the PiezoImplant Passport with them for any restorative procedure. Place a copy of the patient labels from selected abutments in the appropriate location in the PiezoImplant Passport. Please refer to the Rex Implants Restorative Manual for additional restorative instructions.




SP-01 Rev 4, issued January 17, 2019. Additional information is available at <http://www.reximplants.com>.

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Rex Implants, Inc.